



How to set up for iPad or Android Tablet

Once a code has been generated, it will remain inactive until its first time being used, which will then begin the treatment block.

We recommend that the auditory training games are used around 2 times per week for both *Zoo Caper* and *Insane Earplane*. They are meant to last about 15-20 minutes per session. The games will self-timeout at an appropriate time and an internet connection must be used.

Go to either the Apple App Store or Google Play Store and install the correct app. To begin the session, you will need to enter the code generated by the professional into the app. If the professional has transferred the payment of the code to you, then please follow the steps through the PayPal server to complete the activation of the code. Once the code has been purchased, and the app is installed, enter the prescribed code into the app on your iPad or tablet and tap on validate. That will begin the auditory training.

THE TRAINING MUST BE PERFORMED UNDER STEREO HEADPHONES THAT COVERS THE EARS SUFFICIENTLY AND ARE COMFORTABLE TO USE. They do not need to be expensive top-of-the-line models, however ear-buds such as iPod in-ear-buds are not ideal – they let too many low frequencies leak out of the ear. Please set the volume of the iPad or tablet at an easy to hear level that is still quite comfortable for listening.

Your professional will have access to the progress of the auditory training exercise if they want to review the progress.

Once the training is completed, or if there are significant frustrations occurring with the therapy, we recommend you contact the professional for a follow-up evaluation or discussion for any further recommendations.