



Audiometric Testing Suite

Included in the suite is 12 months access to the following:

- New & calibrated audiometric headphones (ANSI/IEC)
- Unlimited air-conduction testing using the Audiometer iPad app which has
 - Automated screening options
 - Manual screening options
 - Both pure-tone and speech screening options (speech testing is more robust in less-than-ideal testing environments)
 - Manual threshold searching capabilities.
- All data is securely sent to our encrypted HIPAA compliant database.
- Unlimited storage of test results in our database which also allows for easy reporting using any specific filtering & tracking with pass/fails of students across any specific time frame.

1-Twelve-month subscription of the Audiometric Testing Suite

€375.00 + shipping

(Cost Does NOT include a required fully functional and supported Apple iPad)



acousticpioneer 

Designed for a superior testing experience



Calibrated Audiometric Headphones

All data is securely stored on HIPAA compliant AWS servers behind multiple encryptions while at rest and in transit. You can easily view results and filter/export data as needed.

acousticpioneer

- Dashboard
- Account
- Users
- Account Credit
- Auditory Therapies
- APD Assessments
- Hearing Screening
- Help

SIGN OUT

Back to home

Hearing Screening

Search...

The search will include the code and patient's name.

Hearing screenings tested areas:

Tones

Speech

Result

Passed

Failed

Threshold

Created at:

From:

to:

[EXPORT TO EXCEL](#) [EXPORT TO CSV](#)

This list shows assessments from the past 30 days. Use the search box to search for specific assessments, including those past 30 days old.

Name	Patient DOB	Created at	Age	User	Right 0.25k	Right 0.5k	Right 1k	Right 2k	Right 3k	Right 4k	Right 6k	Right Speech	Left 0.25k	Left 0.5k	Left 1k	Left 2k	Left 3k	Left 4k	Left 6k	Left Speech	Actions	
emily	10/10/1992	09/28/2022	30	PLEASE	30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass					30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass						
Lisa B	06/30/1987	09/19/2022	35	PLEASE	30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass					30dB - Pass	20dB - Fail	20dB - Fail	20dB - Fail						
Kathleen	03/01/2010	09/12/2022	12	PLEASE	30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass					30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass						
SB	05/03/2013	09/06/2022	9	PLEASE	30dB - Pass	20dB - Fail	20dB - Pass	20dB - Pass					30dB - Fail	20dB - Pass	20dB - Pass	20dB - Pass						
audree	02/04/1961	08/18/2022	61	PLEASE	30dB - Pass	20dB - Pass	20dB - Fail	20dB - Fail					30dB - Pass	20dB - Pass	20dB - Fail	20dB - Fail						
Taylor	08/04/1993	08/11/2022	29	PLEASE	30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass					30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass						
bill	01/10/2007	07/29/2022	15	PLEASE	30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass					30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass						
shoshana	10/12/1984	07/22/2022	37	PLEASE	30dB - Fail	20dB - Pass	20dB - Pass	20dB - Pass					30dB - Fail	20dB - Pass	20dB - Pass	20dB - Pass						
Anna	10/23/1972	07/20/2022	49	PLEASE	30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass					30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass						
Hayley	06/04/1992	07/07/2022	30	PLEASE	30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass					30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass						
old audiologist	12/02/1966	06/07/2022	55	PLEASE	30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass					30dB - Pass	20dB - Pass	20dB - Pass	20dB - Pass				20dB - Fail		

